

Dear Hearing Professional,

I would like to be able to listen in a Hearing Loop. Can you tell me if I have a T-Coil?

- 👂 **If I have a t-coil**, please program it with a Manual T program and/or M/T program set for loop listening. An auto-coil will not work in loop settings so please give me a Manual setting that I can control.
- 👂 **If I do not have a t-coil**, can you tell me if it is possible to have one installed or added with a remote/streamer? If not, please give me options of purchasing a hearing aid with a t-coil.



Assist
A2Hear

Contact Us For Any Questions

Assist2Hear.com
(877) 338-1084